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Instruction Manual

Body Fat Analyzer and Scale

How the Body Fat Analyzer and Scale Works:

This instrument uses BIA (Bioelectrical Impedance Analysis). BIA is one of the most accurate methods of measuring body fat. It simply sends a harmless signal through the fat and muscle in your body. This method calculates both your personal weight and body fat simultaneously, thus giving you a more accurate reading of your overall health and fitness.

Why it is Important to Measure Your Body Fat:

Body fat is a vital measure of your overall health and fitness. Research has shown that healthy body fat levels can help reduce your risk of developing serious health related problems, such as high blood pressure, high cholesterol, heart disease and diabetes. In the past, you would have to consult a medical or fitness professional to obtain a measure of your body fat through an expensive and time-consuming process. Today, with HoMedics Body Fat Analyzer and Scale, simply enter your personal data and step onto the platform.

Caution

BIA (Bioelectrical Impedance Analysis) method determines your body fat percentage by sending a harmless signal through the body. **Do not use this product if you have a pacemaker or other internal medical device. When in doubt, contact your physician.**

This analyzer will give accurate body fat readings for a majority of people, but is not intended for use by the following groups:

Children: Anyone under the age of 10 years

Pregnant Women

Athletes: Anyone that does more than 10 hours or more per week of aerobic activity and has a resting heart rate of 60 beats per minute, especially Body Builders and Professional Athletes.

Note: This instrument is a personal monitor and should be used in a consistent manner for the most accurate reading. While readings might vary between this analyzer and other methods of measurement (hydrodensitometry or hand held calipers), your **changes** in body fat percentage shown by this instrument **will** be reflected accurately.

While this analyzer will work accurately as a scale for anyone and give accurate body fat readings for a majority of people, it is not intended for use by pregnant women, children under 10, or athletic individuals. (see Caution)


Battery Installation

1. Uses one 9-Volt alkaline battery (not included).
2. Remove Battery Compartment Cover on the bottom of the scale and insert battery into the 9-Volt battery snap.
3. The scale needs to be initialized before first use or after battery replacement:
Press lightly on the scale platform. The display will show "0000", then the scale will automatically turn off. Your scale is now ready for use.
4. Replace the battery when "Lo" is displayed.
5. Remove the battery if the product is not used for a long time.

Instructions for Use

Weight Reading Only

1. Select (lb/in) or (kg/cm) by the slide switch on the bottom of the scale.
2. Place the Analyzer/Scale on a firm, flat surface to achieve best accuracy.
3. Step onto the scale; stand still while the scale computes your weight.

4. Your weight will be shown on the LCD display as shown below.

5. The scale turns off by itself after use.


Body Fat Measuring Operation

You must have bare feet to make this measurement. To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on scale. **You must remove your shoes and socks before proceeding with Step 5.** You have approximately 5 seconds to step on the scale in step 5 before it automatically shuts-off to conserve the battery.

1. **Before stepping on scale.** Press (ON) button to enter BFA mode. The height figure of the previous setting will be displayed on LCD. The height figure will be blinking.

2. **Input your height.** Press (▲)/(▼) button to increase /decrease the displayed height. Press (MODE) button to confirm displayed height.

3. **Input your age.** The age figure of the previous setting will be displayed on LCD and it will be blinking.


Press (▲)/(▼) button to increase /decrease the displayed age. Press (MODE) button to confirm displayed age.

4. **Input your gender.** The gender icons will be blinking.

Press (▲)/(▼) button to toggle between male and female icons.


Press (MODE) button to confirm all the settings with 1 long beep.


5. "B" icon will be displayed. Step on the scale with bare feet. Stand still while the unit measures your weight. The weight will display for 2 seconds. Remain on the scale.

6. Then it will start to measure body fat. The measurement pattern will be shown on the display during BFA measurement.

7. The BFA measurement finishes with 1 long beep. The fat percentage, weight, age, and sex will be shown on the display for 4 seconds, then turn off.


- If the product is not used after entering BFA (Body Fat Analyzer) mode, it will automatically turn off in 10 seconds.
- The "Body Fat Display Graph" at the bottom of the screen will show your general condition. See "General Information About Body Fat" for Recommended Body Fat Levels.

Body Fat Measuring Operation Using Memory (optional)

This product features 10 personal memory settings, allowing users to store height, age, and gender for frequent use.

To Store Memory:


- Follow steps 1-5 in Body Fat Measuring Operation.
- Press (▲)/(▼) to select the memory location (0-9) and the memory number will be blinking.

- If the selected memory location is occupied, it will display "■" next to the memory number. If unoccupied, the bar will not be displayed.


- Press (MODE) to store the BFA setting into the desired memory location with 1 long beep.


To Recall Memory:

- Press and hold (ON) button for 2 seconds while the power is off. It will now enter recall mode.
- In Recall mode, it will flash the memory number on LCD for 1 second. The height, age and gender of the memory number will be displayed on the screen.


- Press (▲)/(▼) to select the memory number.
- If the selected memory number is unoccupied, no previous input will be displayed.
- If the selected memory number is occupied, press (MODE) to confirm the desired memory settings with 1 long beep.


- Step on the scale with bare feet. Stand still while the unit measures your weight. The weight will display for 2 seconds.

- Then it starts to measure body fat. The measurement pattern will be shown on the display during BFA measurement.

- The BFA measurement finishes with 1 long beep. The fat percentage, age, and sex will be shown on the display, toggling in 1 second intervals. The unit will turn off automatically after 10 seconds.

To Clear Memory:

- Follow steps 1-4 in **To Recall Memory**.
- If the selected memory is occupied, press both (▲)/(▼) at the same time to clear the content of the selected Memory.

General Information about Body Fat

- Ideal body fat content is not the same for all people. Age, sex, and height are variables in this measurement. The Homedics "Body Fat Display Graph" will show your general condition. Consult your physician to determine what is most ideal for you.

Recommended Body Fat Levels

Age in Years	20-29 yrs.	30-39 yrs.	40-49 yrs.	50-59 yrs.	60+ yrs.
Men	14%	15%	17%	18%	19%
Women	20%	21%	22%	23%	24%

- It is recommended to measure body fat or weight at the same time period daily.
- Use this product without clothing for best accuracy. Clothes can vary in weight and affect the percentage calculation.
- Feet must be bare and clean. For best readings they should also be slightly moist.

Level of Hydration

This analyzer measures your weight, the electrical impedance between your feet, and combines those readings with the information you input (Height, Age, Gender). It then automatically calculates your body fat percentage.

Your level of hydration (water content) varies throughout the day and affects your impedance measurement. If you are dehydrated, there is a greater chance for fluctuation in the body fat percentage displayed by this analyzer than when you are fully hydrated.

The best time for using this analyzer is in the evening. By establishing the pattern of your body's fluctuation over a period of one month, you can determine your average level of body fat. Women will have more fluctuation in water content than men and need to monitor their reading over a longer time period to establish a personal average. You must determine this as a reference point before measuring progress in a health program.

You must wait several hours before taking a body fat reading when: drinking coffee or alcohol, taking diuretic medications, or exercising. These all affect your level of hydration and the accuracy of this analyzer.

Precautions

1. **OVERLOAD WARNING:** remove the weight immediately when the scale displays "Err2"; otherwise, permanent damage to the scale will occur.
2. The product is intended for home/consumer use only; it is not intended for professional use in hospitals or medical facilities.
3. Do not disassemble the product. Other than replacing the battery; it contains no user serviceable parts.
4. Clean after use with a slightly dampened cloth. Do not use solvents or immerse the product in water.
5. Avoid excessive impact or vibration to the unit, such as dropping the product onto the floor.

Trouble Shooting

1. **You must have bare feet when using body fat measurement.** To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on scale.
2. **Display Err0** Initialization error. The user is advised to take the measurement again for an accurate result.
3. **Display Err1** Wrong usage or body not stable. Repeat measurement again.
4. **Display Err2** Overload Warning. Remove the weight immediately; otherwise, permanent damage to the scale will occur.
5. **Display Err3** Body impedance is beyond the technical limit (100 ohm to 1 Kohm). Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.
6. **Display Err4** Body fat percentage is beyond the range (5% to 50%). Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.
7. **Display LB** Replace the battery when "LB" is displayed.
8. The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take the reading with clean, slightly moist feet for best accuracy.

Note: If you are having any problems operating this Body Fat Analyzer please call 800-466-3342, between the hours of 8:30 a.m. to 5:00 p.m. EST or by email at cservice@homedics.com

Specifications

- Uses BIA (Bioelectrical Impedance Analysis) to determine body fat percentage.
- High precision STRAIN GAUGE technology to determine weight.
- Weight Capacity: 330 lb or 150 kg (308 lb or 140 kg for model#SC-550)
- Weight Graduation: 0.2 lb or 100 g
For Model SC-550: 0.2 lb or 100 g up to 220 lb/100 kg
0.5 lb or 200 g over 220 lb
- Body Fat Graduation: 0.1% (0.5% for model#SC-550) (in the range from 4% - 50%)
- Height Range: 3'3.4" - 6'6.7" (100.0 - 199.9 cm)
- Age Range: 10 - 80 years
- Measures both weight and body fat simultaneously in about 10 seconds
- Weight-only feature with automatic step on
- Data can be stored for up to 10 individuals
- User selectable units of measure (lb/in) or (kg/cm)
- Low-battery indicator
- Operates on one 9-Volt alkaline battery (not included)

HoMedics, Inc., guarantees this product free from defects in material and workmanship for the life of the original purchaser from the date of original purchase, except as noted below.

This HoMedics product warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of HoMedics. This warranty is effective only if the product is purchased and operated in the USA. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under warranty. HoMedics shall not be responsible for any type of incidental, consequential or special damages.

To obtain warranty service on your HoMedics product, either hand deliver or mail the unit and your dated sales receipt (as proof of purchase), postpaid, along with check or money order in the amount of \$5.00 payable to HoMedics, Inc. to cover handling.

Upon receipt, HoMedics will repair or replace, as appropriate, your product and return it to you, postpaid. If it is appropriate to replace your product, HoMedics will replace the product with the same product or a comparable product at HoMedics' option. Warranty is solely through HoMedics Service Center. Service of this product by anyone other than HoMedics Service Center voids warranty.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com

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